


Huron Image Manager



Easily store, manage, view and share whole slide images.


Huron Image Manager is a powerful server solution for digital pathology. Lean to deploy and easy to use, the Image Manager completes the digital pathology workflow by connecting users with digital slide images.

Huron Image Manager is comprised of three components:

- an Image Server specifically designed and optimized for digital pathology, supporting image formats unique to the industry
- a Web Server that works with the Image Server for an intuitive, flexible user interface
- an HTML5 Web Viewer

Key Features:

- Zero footprint web viewer – no client software required
- Compatible with all major Windows/Mac browsers
- Annotations including measurements
- Directlink™ image sharing
- Data secured through User Credential Verification
- Supports major file formats in digital pathology



Optimized Viewing and Sharing Experience

- Smoothly view images in-browser using HTML5 viewer
- Navigate using the Slide Thumbnail
- Select the desired magnification with reference to scale
- Annotate and measure extensively
- Instantly share views and collaborate with DirectLink™

